

## Quiltec<sup>®</sup> I Wipe

Two-ply quilted polyester wipe, sealed edge

## Description

The Quiltec I Wipe offers the highest sorptive capacity of any sealed edge wipe available. It is manufactured using two layers of knitted polyester fabric that are laminated together to form a stronger, more adsorbent wipe. Because additional fluid is captured within the quilted pockets between the fabric layers, the Quiltec I Wipe has more than twice the adsorbent capacity of a single layer wipe. Quiltec I has laser-cut, heat-sealed edges and is laundered and packaged within a cleanroom laundry for superior cleanliness characteristics. As a knitted polyester product, it is abrasion and chemical resistant and inherently low in releasable particles and fibers.

Part Numbers are on Page 2

## **Technical Data**

Wipe material	100% knitted polyester			
Wipe construction	2-layer, ultrasonically bonded, Interlock knit			
Attribute; (units)	Value	Test Method		
Basis weight; nominal (g/m²)	240	Contec Method		
Sorbency in water				
Intrinsic; (mL/g)	2.48	IEST-RP-CC004.3, Sec. 8.1		
Extrinsic; (mL/m <sup>2</sup> )	581	IEST-RP-CC004.3, Sec. 8.1		
Sorptive rate; (seconds)	< 1			
Non-volatile residue, NVR		IEST-RP-CC004.3, Sec. 7.1.2		
In deionized water; (g/m <sup>2</sup> )	0.006			
In isopropanol; (g/m²)	0.009			
Specific ions		IEST-RP-CC004.3, Sec. 7.2.2		
Sodium; (ppm)	0.050			
Chloride; (ppm)	0.030			
Particles, readily releasable		IEST-RP-CC004.2, Sec. 5.1		
$P \ge 0.5 \mu m$ ; (x10 <sup>6</sup> /m <sup>2</sup> )	3.25			
Fibers > 100µm; (x10³/m²)	0.220	IEST-RP-CC004.2, Sec. 5.2		

Notes:

1) The information presented here is applicable to the part numbers shown above as well as to any product containing the same materials and produced under the same conditions, regardless of product size or packaging configuration. Please contact a Contec sales representative for more details.

2) Data shown are typical values and should not be used as product specifications.

3) Valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions.

4) Current and/or comparison data may be available. Please contact a Contec sales representative for details.



North America P.O. Box 530 Spartanburg, SC 29304 USA Phone: +1-864-503-8333 www.contecinc.com Europe ZI du Prat – RP 3707 56037 VANNES Cedex FRANCE +33 297 437 690 Revision date: 05/04/15 Page 1 of 2

Copyright © 2015 Contec. Inc. All rights reserved.

NNES Cedex Suzhou 215024 CHINA 137 690 + 86-512-6274 4050

China

No.17 Longyun Road

## Quiltec<sup>®</sup> I Wipe

Part Number	Description	EA/OB1	OB1/OB2	OB2/CS	EA/CS
LWQS0001	Quiltec I Wipe, 9" x 9" (23 x 23cm), sterile	25	4	10	1000
LWQS0002	Quiltec I Wipe, 12" x 12" (30 x 30cm), double bagged, sterile	25	4	5	500
LWQT0019	Quiltec I Wipe, 18" x 18" (46 x 46cm)	25	1	8	200
LWQT0031	Quiltec I Wipe, 12" x 12" (30 x 30cm), sterile	10	5	10	500
LWQT0034	Quiltec I Wipe, 29.5" x 29.5" (75 x 75cm)	20	1	4	80
LWQT0100	Quiltec I Mop Paddle Cover, 10" x 17" (25 x 43cm)	20	1	12	240
QT1-44B	Quiltec I Wipe, 4" x 4" (10 x 10cm), bulk	300	1	8	2400
QT1-66B	Quiltec I Wipe, 6" x 6" (15 x 15cm), bulk	200	1	10	2000
QT1-99	Quiltec I Wipe, 9" x 9" (23 x 23cm)	50	2	10	1000
QT1-99B	Quiltec I Wipe, 9" x 9" (23 x 23cm), bulk	150	1	6	900
QT1-1212	Quiltec I Wipe, 12" x 12" (30 x 30cm)	50	2	5	500
QT1-1212B/100	Quiltec I Wipe, 12" x 12" (30 x 30cm), bulk	100	1	5	500

Part Number	Description	EA/OB1	OB1/OB2	OB2/OB3	OB3/CS	EA/CS
LWQS0004	Quiltec I Wipe, 12" x 12" (30 x 30cm), triple bagged, sterile	25	1	1	20	500

Key: EA = Each, OB1 = Outer Bag 1, OB2 = Outer Bag 2, OB3 = Outer Bag 3, CS = Case

Notes:

1) The information presented here is applicable to the part numbers shown above as well as to any product containing the same materials and produced under the same conditions, regardless of product size or packaging configuration. Please contact a Contec sales representative for more details.

2) Data shown are typical values and should not be used as product specifications.

3) Valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions.

4) Current and/or comparison data may be available. Please contact a Contec sales representative for details.



North America P.O. Box 530 Spartanburg, SC 29304 USA Phone: +1-864-503-8333 www.contecinc.com Europe ZI du Prat – RP 3707 56037 VANNES Cedex FRANCE +33 297 437 690 Page 1 of 2 China No.17 Longyun Road

Suzhou 215024

+ 86-512-6274 4050

CHINA

Revision da

Page 2 of 2